



Quick-Start Operating Guide Document No. 1800-58
Full-Vacuum Control Module

© Copyright 2010 Terra Universal Inc. All rights reserved. • Revised October 2010





Quick-Start Operating Guide Full-Vacuum Control Module

© Copyright 2010 Terra Universal Inc. All rights reserved. • Revised October 2010 • Document No. 1800-58 + 20m

Proprietary Notice

This manual pertains to proprietary devices manufactured by Terra Universal, Inc. Neither this document nor any portion of it may be reproduced in any way without prior written permission from Terra Universal.

Safety Notice

A thorough familiarity with all operating guidelines is essential to safe operation of the product. Failure to observe safety precautions could result in poor performance, damage to the system or other property, or serious bodily injury or death.

The following symbols are intended to call your attention to two levels of hazard involved in operation:

The information presented here is subject to change without notice.

1.0 Introduction

This manual provides information on installing and operating your Terra Full-Vacuum Control Module.

By studying this document carefully, you can be assured of a long, efficient service life from your system.

2.0 Description

Terra Universal's Vacuum Chamber and Vacuum Control Module allow controlled vacuum processing down to one atmosphere (29.9" Hg) before parts are passed into a stainless steel glove box. The turnkey system includes a Vacuum Antechamber with sliding tray (to facilitate parts transport to and from the glove box), a vacuum control module (with vacuum gauge, set point controller and solenoid), and a vacuum pump. Please refer to TUI Operating Manual No. 1800-13 for guidelines on connecting the stainless steel glove box. This QuickStart guide provides instructions on connecting and operating the vacuum control system.



Note: This unit is intended for use with chambers that support a vacuum to 29.9" Hg. For low-vacuum (negative-pressure) chambers, use Terra's Low-Vacuum Controller (No. 1685-99).

3.0 Installation

Carefully unpack all system components and check for missing or damaged parts. Any damage should be reported immediately to the shipping company.

The Vacuum Antechamber is installed on the specified side of the Stainless Steel Glove Box before shipping. Refer to the illustration below for guidelines on connecting the pump and vacuum controller to the antechamber.

Terra Universal makes no warranties applying to information contained in this manual or its suitability for any implied or inferred purpose. Terra Universal shall not be held liable for any errors this manual contains or for any damages that result from its use.



Cautions are used when failure to observe instructions could result in significant damage to equipment.



Warnings are used when failure to observe instructions or precautions could result in injury or death.



4.0 Operation

- 1. Place items to be introduced into the glove box into the vacuum chamber and close both the inside and outside chamber doors.
- 2. Rotate the backfill valve to the OFF position.
- 3. Plug the vacuum pump power cord to the "vacuum pump receptacle" on the rear of the controller. Plug the cable from the solenoid valve to the "solenoid receptacle" on the rear of the controller. Plug the "sensor tube" cable on the back of the controller to the pressure sensor on the vacuum chamber. With the power cord of the controller plugged into a 120VAC/60Hz (for No. 1681-75) or 220VAC/50Hz (for No. 1681-76) receptacle, turn both the solenoid and vacuum pump switches off and the main power on.



- 4. Rotate the chamber value to the OPEN position and turn ON the vacuum pump switch on the vacuum controller. At this point, the vacuum pump is on and ready to evacuate the chamber.
- 5. To set the desired vacuum level you must set the limit switch set point on the vacuum gauge that is integrated into the controller. The controller uses the "LO" set point contacts.

On the vacuum gauge, loosen the thumbscrew at the bottom of the front panel and raise the protection cover. The set point adjusting potentiometers are located to the left of the thumbscrew and are labeled "HI" and "LO." Rotate the front panel switch to "LO" and turn the low potentiometer (clockwise to increase the set point) until the desired set point pressure is displayed

6. To start the evacuation process, turn on the vacuum switch on the controller. This opens the solenoid valve and starts reducing the pressure in the chamber. Evacuation continues until the desired vacuum level is reached as determined by the setting made in step 5 above.





7. As the chamber evacuates, the gauge will read the chamber's vacuum in torr (750 torr is the starting atmospheric pressure). When the set point is reached, the controller will de-energize the solenoid and maintain the desired vacuum inside the chamber.

To retrieve the items from the chamber:

- 1. Turn off the solenoid switch; then turn OFF the "VAC pump" switch on the controller.
- 2. Rotate the "backfill valve" on the chamber to allow the environment inside the glove box to relieve the vacuum in the chamber.
- 3. When the gauge reads 750 torr, you may open either of the vacuum chamber doors.



Quick-Start Operating Guide Full-Vacuum Control Module

© Copyright 2010 Terra Universal Inc. All rights reserved. • Revised October 2010 • Document No. 1800-58 v20mm

5.0 Warranty

Products Manufactured by Terra: Terra Universal, Inc., warrants products that it manufactures to be free from defects for a period of 12 months for parts and 90 days for labor, commencing from the date of shipment. Terra's sole responsibility is to repair or replace, at its option, any part of the product that proves defective or malfunctioning during this time limit. In some cases, components incorporated in Terra Universal products are covered by additional warranties from component manufacturers; obtain specific information from Terra sales representatives. This warranty is void if the equipment is abused or modified by the customer, is operated outside Terra's operating instructions or specifications, or is used in any application other than that for which it is specified. This warranty does not include routine maintenance or service procedures, breakage of quartz baths after 60 days, shipping damage, nor damage from misuse, intentional or unintentional abuse, neglect, natural disasters, or acts of God.

Products Manufactured by Others: Terra Universal, Inc., warrants that, to the best of its ability, Terra's representations of products that are manufactured by others reflect the manufacturer's representations, subject to change without notice. Sole warranty for these products is the original manufacturer's warranty that is passed forward to the purchaser and constitutes the customer's sole remedy for these products. Detailed warranties for distributed products are available through Terra sales representatives.

Freight Shortage or Damage: Upon receipt of any equipment from Terra Universal, Inc., customer shall immediately unpack and inspect for damage or shortage. The customer shall not accept a damaged package or a short shipment until the carrier makes a "damage or shortage" notation on both the carrier's and customer's copy of the freight bill or delivery receipt. Service title passes when the shipment is loaded, so customer is responsible for filing and collecting a freight claim. Any replacement products must be ordered and paid for separately. For Terra's "Policy and Procedures for Returning Goods," see Terra's Internet site: www.TerraUniversal.com.

Generally, customers can improve the chance of collecting on a freight claim by following these procedures: 1) formally requesting that the carrier inspect the shipment immediately upon suspecting damage or shortage to verify condition; 2) notifying the carrier upon discovery of concealed damage and requesting an inspection within 15 days of receipt, both in person or phone and following up via mail; 3) keeping the shipment as intact as possible, including retaining original packaging materials and keeping the product as close to the original receiving location as possible; 4) holding salvage for disposition by the carrier.

All Claims: Terra Universal expressly disclaims all other warranties, expressed or implied or implied by statute, including the warranties of merchantability or fitness for intended use. Terra Universal is not responsible for consequential or incidental damages arising out of the purchase or use of the products supplied by Terra Universal. Terra Universal is not liable for damage to facilities, other equipment, products, property or personnel of others, or of their agents, suppliers, or affiliated parties, which is caused or alleged to have been caused by products supplied by Terra Universal. In any event or series of events, Terra Universal's total liability for any and all damages whatsoever is limited to the lesser of the actual damages or the original invoice cost of the items alleged to have caused the damage. The customer's sole and exclusive remedy for any cause of action whatsoever is repair or replacement of the non-conforming products or refund of the actual purchase price, at the sole option of Terra Universal. All claims must be made in writing within 90 days of the date the product was shipped. Any claims not made within this time limit shall be deemed waived by the customer. Terra Universal is not responsible for any additional costs of repair caused by poor packaging or in-shipment damage during return.

Warranty Returns: All warranty returns must be authorized in advance by Terra Universal and approved under an RMA. Unless approved in advance for good reason, all returns must be in original condition, including all manuals, and must be packaged in original packaging materials. All returned goods are to be shipped to Terra Universal, freight prepaid at customer's expense. See Terra's "Policy and Procedure for Returned Goods."

Thank you for ordering from Terra Universal!