

## Chemical Compatibility Charts

These charts are intended as a general guide for various materials and chemicals. They show some of the materials used in Terra's products and chemicals likely to be used with them. Testing is strongly recommended for extreme conditions of use, such as prolonged exposure or immersion, high temperatures and high concentrations. The acids, caustics and salts in this chart are assumed to be in solution. Materials may react differently to the pure substances (glacial acetic acid, for example).

<b>Hazards</b> (Only the primary ones are shown. For example, chlorine is not shown as an asphyxiant because its toxicity will kill you first).	
A	Asphyxiant (gases and vapors only)
C	Corrosive
F	Flammable
O	Oxidizer
T	Toxic
○	No noticeable effect.
◐	Minor effect or slight change in appearance or properties. Test before repeated exposure.
★	No noticeable effect at low concentration and room temperature. Moderate to severe effect at high concentration and/or high temperature. Test before using.
●	Severe effect or degradation, exposure not recommended.

		METALS						
		HAZARDS-	Aluminum	Brass	Bronze	Copper	304 Stainless Steel	316 Stainless Steel
ACIDS	Acetic	C	◐	●	◐	◐	●	★
	Aqua Regia	C	●	●	●	●	●	●
	Chromic	C	●	●	●	●	◐	◐
	Hydrochloric	C	●	●	●	●	●	●
	Hydrofluoric	C	●		◐	◐	◐	◐
	Nitric	C	●	●	○	●	○	○
	Phosphoric	C	◐	●	◐	●	●	●
Sulfuric	C	●		◐	●	●	●	
CAUSTICS	Ammonium Hydroxide	C	◐	●	●	●	○	○
	Potassium Hydroxide	C	●	●	●	◐	◐	○
	Sodium Hydroxide	C	●	●	◐	●	◐	◐
GASES	Air	O	○	○	○	○	○	○
	Ammonia	C,F,T	○	●	●	●	○	○
	Argon	A	○	○	○	○	○	○
	Carbon Dioxide	A	○	○	○	○	○	○
	Carbon Monoxide	F,T	○	○	○	○	○	○
	Chlorine	C,T	●	●	●	●	○	○
	Flourine	C,O,T	★	★	◐	★	○	○
	Helium	A	○	○	○	○	○	○
	Hydrogen	A,F	○	○	○	○	○	○
	Hydrogen Sulfide	C,F,T	○	○			○	○
	Methane	A,F	○	○	○	○	○	○
	Nitrogen	A	○	○	○	○	○	○
	Nitrous Oxide	O	○	○	○	○	○	○
	Ozone	O	◐		◐	○	◐	○
Propane	A,F	○	○	○	○	○	○	
OXIDIZERS	Hydrogen Peroxide	O	○	●	◐	●	◐	○
	Sodium Hypochlorite	O	●	●	◐		★	★
SALTS	Ammonium Nitrate		◐	●	●	●	○	○
	Ammonium Persulfate		●	●	●	●	●	◐
	Sodium Carbonate		●	◐	○	○	○	○
SOLVENTS	Acetone	F	○	○	○	○	○	○
	Carbon Tetrachloride	T	●	○	○		◐	◐
	DI Water		○	○		◐	○	○
	Ethyl Alcohol	F	◐	○	○	○	○	○
	Ethylene Glycol		○	◐	○	○	◐	◐
	Glycerine		○	◐	○	○	○	○
	Isopropyl Alcohol	F	◐		○	◐	◐	◐
	Kerosene	F	○	○	○	◐	○	○
	Methyl Alcohol	F,T	○	○	○	◐	○	○
	Methyl Ethyl Ketone	F	◐	○	○	○	○	◐
Toluene	F	○	○	○	○	○	○	
Trichloroethane	A	●		○		◐	◐	

		PLASTICS																
		-HAZARDS-	Acrylic (plexiglass) ABS	CPVC	Noryl	Nylon	Polycarbonate	Polyethylene	Polypropylene (TPX)	Polypropylene	Polystyrene	PPS (Ryton)	PVC, Static Dissipative PVC	Styrene Acrylonitrile (SAN) PVDF (Kynar)	Tellon, PTFE			
ACIDS	Acetic	C	●	*	*	○	●	*	●	*	*	○	*	*	○	*	○	
	Aqua Regia	C	●		▶	●	●	●	●	●	●	●	●	●	○	*	○	
	Chromic	C	▶	▶	*	*	●	*	*	*	*	●	○	*	○	*	○	
	Hydrochloric	C	*	○	*	○	●	*	*	*	*	*	●	*	○	*	○	
	Hydrofluoric	C	●	▶	●	*	●	●	*	*	○	*	○	*	*	*	○	
	Nitric	C	●	*	*	*	●	*	*	*	*	●	●	*	○	*	▶	○
	Phosphoric	C	●		*	○	●	○	*	○	▶	○	●	▶	○	*	○	
	Sulfuric	C	*	*	*	○	●	*	*	○	*	*	*	*	*	*	○	○
CAUSTICS	Ammonium Hydroxide	C	▶	○	*	○	○	●	*	○	○		○	*	○	*	○	
	Potassium Hydroxide	C	▶		○	○	▶	●	○				*	*	○	*	○	
	Sodium Hydroxide	C	▶	○	*	○	*	*	*	○	▶	○	▶	*	○	*	○	
GASES	Air	O	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	
	Ammonia	C,F,T	▶		*	▶	○	●	●	○	▶		*	○	○	*	○	
	Argon	A	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	
	Carbon Dioxide	A	▶		*	○	○	○	*	○	▶	○	○	*	○	*	○	
	Carbon Monoxide	F,T			*	○	○	○	*		▶		○	*	○	*	○	
	Chlorine	C,T			●	●	●	○	●	●	●		●	●	●	*	○	
	Flourine	C,O,T	○		*		●	▶	●	●	●		○	●	●	*	●	
	Helium	A	○	○	*	○	○	○	○	○	▶	○	○	*	○	*	○	
	Hydrogen	A,F			*	○	○	○	▶		▶		*	*	○	▶	○	
	Hydrogen Sulfide	C,F,T	▶		*		●	▶	*		▶		▶	*		*	○	
	Methane	A,F			*	●	○	○	○		▶		○	*		*	○	
	Nitrogen	A	○	○	○	○	○	○	○	○	○		○	○	○	●	○	
	Nitrous Oxide	O			*		▶		○		●		○	*		○	○	
	Ozone	O	○		○	○	●	●	○		▶			▶		○	○	
	Propane	A,F			○	○	○	▶	○		▶		○	*	○	*	○	
OXIDANTS	Hydrogen Peroxide	O	○	○	*	○	●	○	○	*	*	○	●	●	*	*	○	
	Sodium Hypochlorite	O	▶	○	*	○	●	▶	*	○	*	○	●	*	○	*	○	
SALTS	Ammonium Nitrate		○		*	○	●	○	*	○	*	○	○	*	○	*	○	
	Ammonium Persulfate		○		*	○	●		*		*		○	*	○	*	○	
	Sodium Carbonate		*	○	*	○	○	○	○	○	*	○	*	*	○	*	○	
SOLVENTS	Acetone	F	●	●	●	●	○	●	●		○	●	*	●	●	●	○	
	Carbon Tetrachloride	T	●	▶	●	●	●	●	●	●	●	●	○	●	●	*	○	
	DI Water		○	○	○	○	○	○	○		○		○	○	○	○	○	
	Ethyl Alcohol	F	○	▶	▶	○	○	▶	*		*		○	▶	○	○	▶	
	Ethylene Glycol		▶		*	○	○	○	▶	○	▶	○	*	*	○	*	○	
	Glycerine		▶		*	○	○	○	▶		▶	○	*	*	○	*	○	
	Isopropyl Alcohol	F		▶	▶	○	●	○	*		*			▶	○	*	○	
	Kerosene	F	▶	○	*	●	○	*	*		●		*	▶	○	*	○	
	Methyl Alcohol	F,T	●	▶	*	○	○	▶	*		*		○	▶	○	*	●	
	Methyl Ethyl Ketone	F	●		●	●	▶	●	●	▶	*	●	*	●	●	●	●	○
	Toluene	F	●	●	●	●	○	●	●	●	●	●	*	●	●	*	●	○
Trichloroethane	A				●	▶	●	○		▶			●	○	○	○	○	

RUBBER & SYNTHETICS (ELASTOMERS)												
-HAZARDS-												
Buna N												
EPDM												
Hypalon												
Natural Rubber												
Neoprene												
Nitrile												
PVA (Polyvinyl Alcohol)												
Tygon												
Silicone												
Viton												
ACIDS	Acetic	C	▶	○	*	*	*	▶	●	●	▶	*
	Aqua Regia	C	●	▶	▶	●	*	*	●	●	●	*
	Chromic	C	●	▶	▶	●	●	▶	●	▶	▶	*
	Hydrochloric	C	▶	*	*	○	*	▶	●	●	●	*
	Hydrofluoric	C	●	●	▶	*	●	○	●	●	●	*
	Nitric	C	●	*	*	●	●	*	●	●	●	*
	Phosphoric	C	●	▶	▶	▶	*	*	●	●	●	*
	Sulfuric	C	*	*	*	*	*	●	●	●	●	*
CAUSTICS	Ammonium Hydroxide	C	○	○	○	*	*	○	●	○	○	●
	Potassium Hydroxide	C	▶	○	○	○	○	○	●	▶	▶	▶
	Sodium Hydroxide	C	*	▶	○	○	*	○	●	▶	○	*
GASES	Air	O	○	○	○	○	○	○	○	○	○	○
	Ammonia	C,F,T	▶	○	●	●	○	○		○	▶	●
	Argon	A	○	○	○	○	○	○	○	○	○	○
	Carbon Dioxide	A	*	*	*	*	*	*			*	*
	Carbon Monoxide	F,T	○	○	▶	●	*	*			○	*
	Chlorine	C,T	▶	○	●	●	●			▶	●	*
	Flourine	C,O,T	●	○		▶	●				●	*
	Helium	A	○	○	○	○	*	○	○	○	○	○
	Hydrogen	A,F	○		○	*	*			○	▶	*
	Hydrogen Sulfide	C,F,T	●	▶	▶	▶	●			●	▶	*
	Methane	A,F	○	●	▶	●	*				●	*
	Nitrogen	A	○	○	○	○	*	○		○	○	*
	Nitrous Oxide	O	○				●					*
	Ozone	O	●	○	○	●	▶				○	*
	Propane	A,F	○	●		●	*				●	*
OXIDANTS	Hydrogen Peroxide	O	●	*	●	▶	●	○	●	▶	*	*
	Sodium Hypochlorite	O	▶	▶	*	▶	*		●	▶	▶	*
SALTS	Ammonium Nitrate		○	○	○	▶	▶	○	●	○	▶	*
	Ammonium Persulfate		○	▶	○	○	*	○	●	○	●	*
	Sodium Carbonate		○	○	○	○	*		●	▶	○	*
SOLVENTS	Acetone	F	●	○	▶	▶	*	*	●	●	▶	●
	Carbon Tetrachloride	T	▶	▶	●	●	*	*	○		●	*
	DI Water		○	○	○	○	○		●	○	○	○
	Ethyl Alcohol	F	▶	○	○	○	*	○	●	▶	▶	*
	Ethylene Glycol		○	○	○	○	*	○	▶	▶	○	*
	Glycerine		○	○	○	○	*	○	▶	○	○	*
	Isopropyl Alcohol	F	▶	○	○	○	*	○	●	○	○	*
	Kerosene	F	○	●	●	●	●	○	▶	●	●	*
	Methyl Alcohol	F,T	▶	○	○	○	*	○	●	○	▶	*
	Methyl Ethyl Ketone	F	●	○	●	*	●	●	▶	●	●	●
	Toluene	F	●	●	●	●	*	*	▶	●	●	*
Trichloroethane	A	●	●	●	●	*	*	▶		●	○	