

## Chemical Filtration Table

The table below indicates the filter recommended for use with a wide range of chemicals commonly found in laboratory and manufacturing environments. For each chemical, it also provides an estimate of filter capacity.

Unless otherwise specified, filter efficiency with the chemical indicated approaches 100%. An asterisk (\*) indicates poor adsorption by any filter; use of the specified filter is not recommended with these chemicals except where they occur in low concentrations, and frequent exhaust testing is advised.

Filter capacity is given as the equilibrium saturation capacity, a standard test procedure for activated carbon (ASTM-D3467). The capacity is expressed as the final filter weight at saturation as a percentage of initial filter weight. The actual capacity obtained will depend on the conditions of use (concentration of contaminant feed, temperature, humidity, etc.).

**Note:** All data below represent typical, not guaranteed, test results. Terra recommends thorough testing of filter performance in each application prior to exhaust release, and assumes no responsibility for damage or injury that may result as a consequence of improper filter selection or use.

| Chemical                      | Filter | Satur. Cap. (%) |
|-------------------------------|--------|-----------------|
| <b>Acids</b>                  |        |                 |
| Acetic                        | AC     | 33              |
| Acetic anhydride              | AC     | 33              |
| Acrylic                       | AC     | 40              |
| Butyric                       | AC     | 40              |
| Caprylic                      | AC     | 40              |
| Carbolic                      | AC     | 40              |
| Formic                        | AC/I   | 20              |
| Lactic                        | AC     | 40              |
| Osmium tetroxide              | AC     | 40              |
| Palmitic                      | AC     | 40              |
| Phenol                        | AC     | 40              |
| Propionic                     | AC     | 40              |
| Valeric                       | AC     | 40              |
| Alcohols                      | x      | x               |
| Amyl                          | AC     | 40              |
| Butyl                         | AC     | 40              |
| Cyclohexanol                  | AC     | 45              |
| Ethyl                         | AC     | 32              |
| Isopropyl                     | AC     | 40              |
| Methyl (methanol)             | AC     | 32              |
| Propyl                        | AC     | 40              |
| <b>Aliphatic Hydrocarbons</b> |        |                 |

| Chemical                       | Filter | Satur. Cap. (%) |
|--------------------------------|--------|-----------------|
| <b>Aromatic Hydrocarbons</b>   |        |                 |
| Benzene                        | AC     | 40              |
| Napthalene                     | AC     | 47              |
| Ninhydrin                      | AC     | 47              |
| Styrene monomer                | AC     | 47              |
| Toluene                        | AC     | 47              |
| Toluidine                      | AC     | 47              |
| Xylene                         | AC     | 40              |
| <b>Esters</b>                  |        |                 |
| Butyl acetate                  | AC     | 40              |
| Cellosolve acetate             | AC     | 45              |
| Ethyl acetate                  | AC     | 40              |
| Ethyl acrylate                 | AC     | 45              |
| Ethyl formate                  | AC     | 40              |
| Isopropyl acetate              | AC     | 45              |
| Methyl acetate                 | AC     | 40              |
| Methyl acrylate                | AC     | 45              |
| Methyl formate                 | AC     | 40              |
| <b>Aldehydes &amp; Ketones</b> |        |                 |
| Acetone                        | AC     | 32              |
| Acetaldehyde                   | FOR    | 10              |
| Acrolein                       | AC     | 32              |
| Benzaldehyde                   | AC     | 40              |

|             |    |    |
|-------------|----|----|
| Acetylene   | AC | 20 |
| Iso-butane  | AC | 10 |
| Butylene    | AC | 10 |
| Butadiene*  | AC | -- |
| Cyclohexane | AC | 35 |
| N-decane*   | AC | -- |
| Ethane*     | AC | -- |
| Ethylene*   | AC | -- |
| N-heptane*  | AC | -- |
| Heptylene*  | AC | -- |
| Hexane      | AC | 35 |
| Hexylene*   | AC | -- |
| Methane*    | AC | -- |
| N-nonane*   | AC | -- |
| N-octane*   | AC | -- |
| N-octylene* | AC | -- |
| Pentane     | AC | 26 |
| Propane*    | AC | -- |
| Propylene   | AC | 10 |

|                        |      |    |
|------------------------|------|----|
| Butyraldehyde          | AC   | 32 |
| Caproaldehyde          | AC   | 40 |
| Crotonaldehyde         | AC   | 40 |
| Cyclohexanol           | AC   | 40 |
| Diethyl ketone         | AC   | 32 |
| Dipropyl ketone        | AC   | 40 |
| Formaldehyde           | FOR  | 10 |
| Mesityl oxide          | AC   | 40 |
| Methyl butyl ketone    | AC   | 40 |
| Methyl ethyl ketone    | AC/I | 32 |
| Methyl isobutyl ketone | AC   | 40 |
| Propionaldehyde        | AC   | 32 |
| Valeraldehyde          | AC   | 40 |
| Valeric aldehyde       | AC   | 40 |
| <b>Ethers</b>          |      |    |
| Amyl                   | AC   | 35 |
| Butyl                  | AC   | 35 |
| Cellosolve             | AC   | 40 |
| Dioxan                 | AC   | 45 |
| Diethyl (ethyl)        | ETH  | 10 |
| Ethylene oxide         | AC   | 20 |
| Isopropyl              | AC   | 25 |
| Methyl                 | ETH  | 10 |
| Methyl cellosolve      | AC   | 45 |
| Propyl                 | AC   | 30 |

| Chemical                | Filter | Satur. Cap. (%) |
|-------------------------|--------|-----------------|
| <b>Halogens</b>         |        |                 |
| Bromine                 | AC     | 53              |
| Butyl chloride          | AC     | 40              |
| Carbon tetrachloride    | AC     | 65              |
| Chlorine                | AC     | 20              |
| Chlorobenzene           | AC     | 53              |
| Chlorobutadiene         | AC     | 40              |
| Chloroform              | AC     | 60              |
| Chloro picrin           | AC     | 65              |
| Chloro nitropropane     | AC     | 60              |
| Dibromoethane           | AC     | 60              |
| Dichlorobenzene         | AC     | 60              |
| Dichlorodifluoromethane | AC     | 20              |
| Dichlorodifluoroethane  | AC     | 40              |
| Dichloro ethyl ether    | AC     | 53              |
| Dichloromethane         | AC     | 53              |
| Dichloromono-           | AC     | 20              |

| Chemical                  | Filter | Satur. Cap. (%) |
|---------------------------|--------|-----------------|
| <b>Nitrogen Compounds</b> |        |                 |
| Ammonia                   | AM     | 10              |
| Amines low MW             | AM     | 10              |
| Amines high MW            | AM     | 40              |
| Aniline                   | AC     | 40              |
| Diethyl amine             | AM     | 20              |
| Diethyl aniline           | AC     | 53              |
| Dimethyl amine            | AM     | 20              |
| Ethyl amine               | AM     | 20              |
| Hydrogen cyanide          | CYN    | 20              |
| Indole                    | AC     | 53              |
| Nicotine                  | AC     | 40              |
| Nitric acid fumes         | AC/I   | 10              |
| Nitrobenzene              | AC     | 53              |

|                            |      |    |
|----------------------------|------|----|
| fluoromethane              |      |    |
| Dichloropropane            | AC   | 53 |
| Dichlorotetra-fluoroethane | AC   | 20 |
| Ethyl bromide              | AC   | 20 |
| Ethyl chloride             | AC   | 20 |
| Ethylenechlorohydrin       | AC   | 40 |
| Ethylene dichloride        | AC   | 53 |
| Fluorotrichloro-methane    | AC   | 50 |
| FREON®(BP> -20°C)          | AC   | 45 |
| Hydrogen bromide           | AC/I | 5  |
| Hydrogen chloride          | AC/I | 4  |
| Hydrogen iodide            | AC/I | 7  |
| Iodine                     | AC   | 55 |
| Iodoform                   | AC   | 53 |
| Methyl bromide             | AC   | 25 |
| Methyl chloride            | AC   | 20 |
| Methyl chloroform          | AC   | 45 |
| Methylene chloride         | AC   | 45 |
| Monochlorobenzene          | AC   | 45 |
| Paradichlorobenzene        | AC   | 45 |
| Perchloroethylene          | AC   | 45 |
| Propyl chloride            | AC   | 40 |
| Tetrachloroethane          | AC   | 53 |
| Tetrachloroethylene        | AC   | 53 |
| Vinyl chloride             | AC   | 20 |
| <b>Sulfur Compounds</b>    |      |    |
| Carbon disulfide           | AC   | 20 |
| Dimethyl sulfate           | AC   | 50 |
| Ethyl mercaptan            | SUL  | 40 |
| Hydrogen sulfide           | SUL  | 20 |
| Mercaptans high MW         | SUL  | 40 |
| Sulfur dioxide             | AC/I | 10 |
| Sulfur trioxide            | AC/I | 20 |
| Sulfuric acid              | AC/I | 40 |
| Tetrahydrothiapene         | AC   | 40 |

|                      |      |    |
|----------------------|------|----|
| Nitroethane          | AC   | 53 |
| Nitrogen dioxide*    | AC/I | —  |
| Nitromethane         | AC   | 40 |
| Nitropropane         | AC   | 40 |
| Nitrotoluene         | AC   | 53 |
| Pyridine             | AM   | 53 |
| Urea                 | AC   | 53 |
| Uric acid            | AC   | 53 |
| <b>Miscellaneous</b> |      |    |
| Adhesives            | AC   | 40 |
| Animal Odors         | ALK  | 30 |
| Camphor              | AC   | 40 |
| Carbon monoxide*     | AC   | —  |
| Carbon dioxide*      | AC   | —  |
| Citrus fruits        | AC   | 40 |
| Cooking odors        | AC   | 40 |
| Deodorizers          | AC   | 20 |
| Detergents           | AC   | 40 |
| Hospital odors       | ACD  | 30 |
| Human odors          | ACD  | 30 |
| Leather              | AC   | 30 |
| Ozone                | AC   | 30 |
| Perfumes             | AC   | 30 |
| Gasoline             | AC   | 40 |
| Putrescine           | ACD  | 30 |
| Resins               | AC   | 30 |
| Toilet odors         | ALK  | 30 |

| <b>Filter Types</b> |                             |     |         |
|---------------------|-----------------------------|-----|---------|
| ACD                 | Organic Acid                | SUL | Sulphur |
| AC/I                | Inorganic Acid              | ETH | Ether   |
| ALK                 | Alkaline Odor               | CYN | Cyanide |
| FOR                 | Formaldehyde                |     |         |
| AC                  | Activated Carbon            |     |         |
| AM                  | Amines (nitrogen compounds) |     |         |