






# XPert<sup>®</sup> Balance and Bulk Powder Enclosures




# Labconco Product Diversity



# Fume Hood or Enclosure?




Chemical Fume Hood



Balance Enclosure

# Fume Hood or Enclosure?

## Chemical Fume Hoods



- Capture, contain and remove airborne contaminants
- Large energy consumption
- Low duct velocity, 2000-2500 fpm
- High installation cost

# Fume Hood or Enclosure?

## Balance Enclosure

- Capture, contain and remove airborne particulate contaminants
  - Small energy consumption
- High duct velocity, 3500-4000 fpm
  - HEPA filter required
  - Low installation cost



# XPert® Balance Enclosures and Systems

## Markets

- Pharmaceutical
- Industrial Hygiene
- Chemistry Labs
  - Life Sciences
  - University
- Hospital/Clinical




# XPert<sup>®</sup> Filtered Balance Systems



- Provide user protection
- Patented rear baffle with zones of perforation
- Built-in blower – need no ducting to outside
- Exhaust HEPA filter
- Isolation supports minimize blower vibration
- Low installation costs
- Available in 2', 3', 4', 5', and 6' widths
- Optional left-side mounted waste chute



# XPert<sup>®</sup> Filtered Balance Systems

## Features & Benefits





# XPert<sup>®</sup> Filtered Balance Systems



- Complete “system”—self contained
  - ✓ Internal fan w/vibration isolation supports
  - ✓ Filter mounted above integral fan
  - ✓ Fits on existing bench
- Exhaust options
  - ✓ Recirculate into room
  - ✓ Exhaust to house via canopy


# XPert® Filtered Balance Stations




- Protect the user
- Built-in HEPA filter removes particulates
- Must be ducted to the outside using a remote blower
- Patented rear baffle
- Powder weighing applications in pharmaceutical research are typical
- Deep interior accommodates large micro and analytical balances
- Optional left-mounted waste chute

# XPert® Filtered Balance Stations

## Features & Benefits



# XPert® Filtered Balance Stations



- No internal fan
- Filter mounted above enclosure
- Tie into house exhaust system
- Minimized vibration

# XPert<sup>®</sup> Balance Enclosures



- Provide user protection
- Must be ducted to the outside or connected to a FilterMate Portable Exhauster
- Deep interior accommodates large micro and analytical balances
- Patented rear baffle
- Exhaust less tempered room air than fume hoods — conserving energy
- No switches, wiring or service fixtures that can add costs


# XPert<sup>®</sup> Balance Enclosures



- 2', 3', 4', 5', 6', & 8' models
- 22.3" or 31.5" exterior height
- 23.4" or 30.0" interior depth
- Top or bottom exhaust connection
- Airflow monitor standard
- Optional left-side mounted waste chute

# XPert<sup>®</sup> Balance Enclosures

## Features & Benefits



# XPert® Balance Enclosures

## Easy-Order




### One part number includes:

- 2', 3' or 4' Balance Enclosure
- FilterMate Portable Exhauster
- Bag-In/Bag-Out HEPA Filter
- Vent Connection Kit with adapters and hose connectors
- Guardian Jr. or 1000 Digital Airflow Monitor
- Dished Epoxy Work Surface



# XPert<sup>®</sup> Balance Enclosures


## Performance-Enhancing Features



- Enhanced Containment
  - ✓ Horizontal slotted baffle
  - ✓ Clean-Sweep™ air foil
  - ✓ Bypass air
  - ✓ Sash foil
  - ✓ Side-entry air foils


# XPert<sup>®</sup> Balance Enclosures

## Containment-Enhancing Features



# XPert<sup>®</sup> Balance Enclosures

## Containment-Enhancing Features



**Clean-Sweep™ Air Foil** allows air to sweep the work surface for maximum containment.



**Side-Entry Air Foils** are located on either side of the sash opening and direct clean airflow into the chamber.



**Upper Dilution Air Supply** introduces air from above the work area and bathes the back of the sash with clean air and directs concentrations away from the sash opening.

# XPert<sup>®</sup> Balance Enclosures

## Containment-Enhancing Features



### Aerodynamic Air Foil

- Clean-Sweep<sup>™</sup> openings
- Directs airflow into hood and across work surface
- Minimizes turbulence
- Patented feature




# XPert<sup>®</sup> Balance Enclosures

## Containment-Enhancing Features

### Upper Containment Sash Foil

- Aerodynamic shape directs airflow
- Air gap bleeds high velocity air into enclosure
- Reduces contaminants behind sash
  - Patented feature




# XPert<sup>®</sup> Balance Enclosures

## Containment-Enhancing Features

### Side-Entry Air Foils

- Provides horizontal flow
- Eliminates dead-air zone around sash area
  - Exclusive feature



# XPert<sup>®</sup> Balance Enclosures

## Containment-Enhancing Features



### Upper Dilution Air Supply

- Bypass Air (5 – 10%)
- Purges top of enclosure
- Bathes area behind sash
  - ✓ Reduces chemical concentrations
  - ✓ Introduces clean air near user's breathing zone
- Exclusive feature

# XPert<sup>®</sup> Balance Enclosures

## Containment-Enhancing Features

### Slotted Horizontal Baffle

- Greatly reduces vortex (roll)
- Provides horizontal “laminar” flow
  - Patented feature





# XPert® Balance Enclosures

## Vibration Isolation



### Vibration:

- Affects weighing operations
- Greater effect in taller buildings
- Placement near outside wall is better
- Corners are more rigid
- Work surface needs to be stable and level

# XPert<sup>®</sup> Balance Enclosures

## Validation Documentation

- Validation package included with each enclosure:
  - ✓ Independent test data:
    - SafeBridge Consultants: Powders Containment Test
    - Exposure Control Technologies: ASHRAE 110-1995 Test
  - ✓ IQ/OQ Protocol
  - ✓ Manufacturer's quality assurance inspection form



# XPert<sup>®</sup> Balance Enclosures

## Factors Affecting Performance


- Electrostatics
- Temperature
  - Humidity
  - Vibration
- Lab air currents
- Exhaust system
  - Installation
- User technique



# XPert<sup>®</sup> Balance Enclosures

## Electrostatics

- Electrostatic Problems:



- ✓ Lightweight powders float

- ✓ Static attraction can be greater than air currents

- ✓ Balance performance can be adversely impacted

- Electrostatic problems are compounded by:

- ✓ Enclosure materials of construction

- ✓ Low humidity—dry air holds static charge

- ✓ Temperature—heat sources affect air currents

# XPert<sup>®</sup> Balance Enclosures

## Electrostatics


- Tempered safety glass is static dissipative – releases static more readily
- Acrylic has insulative properties – static charge clings to material
  - Recommend the use of ionizers – aligns charges



Material	Surface Resistivity	Example
Conductive	$0 \rightarrow 10^5 \Omega$ per square CM	Skin, Metals
Static dissipative	$10^5 \rightarrow 10^9 \Omega$ per square CM	Glass
Antistatic	$10^9 \rightarrow 10^{12} \Omega$ per square CM	Polyethylene bag
Insulative	$10^{12} \Omega$ per square CM $\rightarrow$	Acrylic box, Packing foam, Styrofoam

# XPert® Balance Enclosures

## Vibration Isolation




Marble Slab

- 3" Marble Slab
  - ✓ Economical
  - ✓ Fits comfortably
  - ✓ Does not affect containment
- Rubber/Vinyl Isolator Pads
  - ✓ Inexpensive
  - ✓ Effective vibration dampening

# XPert<sup>®</sup> Balance Enclosures


## HVAC/Installation



- Installation Location:
  - ✓ Low traffic areas
  - ✓ Away from open doors
  - ✓ No cross-drafts
  - ✓ Avoid supply air diffuser
  - ✓ Lab location
  - ✓ Low “throw velocity”
- Balanced HVAC system
  - ✓ Adequate supply air

# XPert® Balance Enclosures


Typical venting with FilterMate  
Portable Exhauster






# XPert<sup>®</sup> Balance Enclosures

Typical venting from top with remote blower



# XPert® Balance Enclosures

Typical venting from back with remote blower



# XPert® Balance Enclosures

## FilterMate™ Portable Exhauster

- Point-of-use exhauster
- Internal spark-proof impeller
  - Houses filter to purify air
  - Provides airflow to XPert



# XPert<sup>®</sup> Balance Enclosures

## FilterMate<sup>™</sup> Portable Exhauster

- Narrow design fits comfortably on bench or on the floor
  - Comes standard with 8' polypropylene hose, clamps and Exhaust Transition Adapter



# XPert<sup>®</sup> Balance Enclosures

## FilterMate<sup>™</sup> Portable Exhauster



- Large capacity – 280 CFM
  - ✓ Provides up to 100 fpm face velocity for 2', 3' or 4' Enclosure\*\*
  - ✓ HEPA Filter, Carbon Filter or BOTH with one FilterMate
  - ✓ Carbon filter types:
    - Organic
    - Ammonia
    - Formaldehyde
- Quiet operation
  - ✓ 54 dbA at 150 CFM operation

\*\* Contact Labconco for technical assistance regarding using multiple FilterMates for 5', 6' and 8' wide enclosures

**LABCONCO**

# XPert® Balance Enclosures

## FilterMate™ Portable Exhauster



- Three ways to exhaust:
  - ✓ Through FilterMate to house exhaust system
  - ✓ Through FilterMate back into the room
  - ✓ Directly from enclosure to house exhaust system



# XPert<sup>®</sup> Balance Enclosures

## FilterMate<sup>™</sup> Portable Exhauster



- Intrinsically safe design
  - ✓ Filter is in a negative pressure plenum
- Bag-In/Bag-Out HEPA filter replacement
  - ✓ Minimizes servicer contact
- HEPA filter can be scanned for leaks from outside FilterMate

# XPert® Balance Enclosures

## Convenience Features

- Adjustable speed control (qualified personnel only)
  - Additional grounded outlet
    - ✓ Can power airflow monitor
    - Circuit breaker – 5 amp





# XPert<sup>®</sup> Balance Enclosures

## User Technique



- Use glass or metal weigh dishes and tools – static dissipative
- Work slowly, methodically
- Be cognizant of surroundings
- Work in a bright environment; rest eyes often

# XPert® Balance Enclosures

## User Technique

- No enclosure can compensate for poor technique.
  - Labconco's Enclosure has been proven to provide a safe work environment even with less proficient technique.

*“The Labconco XPert Balance Enclosure is capable of providing excellent containment where the operator uses proficiency in technique...”*

*“Where operators exercise less than excellent technique and/or have minor spillage good containment performance is achieved.”*

**Edward J. Haggerty, CIH**  
**SafeBridge Consultants, Inc.**




# XPert<sup>®</sup> Balance Enclosures

## Safety Features



- Spring-loaded sash latch
  - ✓ Exclusive feature
  - ✓ Allows for easy loading of equipment
- Airflow monitor standard
  - ✓ Guardian Jr. (low flow)
  - ✓ Guardian 1000 (high/low flow)
  - ✓ Easy to install
  - ✓ Calibrated in house

# XPert® Bulk Powder Transfer System



- Definitions & Terms
- Features
- Planning & Start Up
- Loading
- Technique
- Procedure
- Clean Up

# XPert® Bulk Powder Enclosures

## Definitions



**Source Container:** Bulk container or fiber drum containing powder

**Transfer Cart:** Transport device for source container

**Transfer Plate:** Platform of Transport Cart that connects to the lift

**Electric Lift:** Raises the source container for access through enclosure

**Transfer Cart Latch:** Secures Transport Cart to electric lift

# XPert<sup>®</sup> Bulk Powder Enclosures

## Definitions



**Safety Sleeve:** Bag placed around outside of source container for additional protection


**Retaining Band:** Secures the Safety Sleeve to the Source Container and secures the source container bags to the retaining ring

**Retaining Ring:** Ring on work surface used to secure source container bags in place

**Waste Chute:** Hole on side of enclosure that allows attachment of bag for waste disposal

# XPert® Bulk Powder Enclosures

## Terms



Transfer Cart Latch

Electric Lift

Source Container

Transfer Plate

Transfer Cart

# XPert® Bulk Powder Enclosures

## Markets



Pharmaceutical Research



Chemical Manufacturing




# XPert® Bulk Powder Filtered Systems, Filtered Stations and Enclosures



- Designed to protect user during powder transfer operations
- Filtered Systems and Stations include a built-in HEPA filter to remove particulates in the exhaust air
- Patented rear baffle with zones of perforations promotes horizontal airflow

# XPert® Bulk Powder Filtered Systems, Filtered Stations and Enclosures



# XPert® Bulk Powder Filtered Systems



- Built-in HEPA filter to remove particulates in the exhaust air
- Built-in blower — **no ducting to the outside required**
- Push-button operated electric hydraulic lift mechanism
- Electronic beam automatically stops lift mechanism at optimal access height
- Available in 5- and 6-foot widths

# XPert<sup>®</sup> Bulk Powder Filtered Systems



- Required accessories include
  - ✓ Retaining Ring
  - ✓ Safety Sleeves
  - ✓ Retaining Bands
  - ✓ Access Opening Cover

# XPert® Bulk Powder Filtered Stations



- Have an exhaust collar and **require ducting to house exhaust or a remote blower**
- Available in 5' and 6' widths
- Require Remote Blower
- Ductwork
- Retaining ring and bands
- Access Opening Cover


# XPert<sup>®</sup> Bulk Powder Enclosures



- Requires venting by connection to a remote blower or house exhaust to the outside
- Available in 5', 6', and 8' widths
- Require Retaining Ring and Bands
- Require Safety Sleeves
- Require Access Opening Cover

# XPert® Bulk Powder Enclosures

## Exclusive Features




- Transfer Cart is latched in place to keep it from moving
- Transfer Cart handle is out of the way during the weigh procedure

# XPert® Bulk Powder Enclosures

## Exclusive Features

- Lift is electric and automatically stops at access height.
- Transfer Plate locates source container in correct position for access opening.






# XPert<sup>®</sup> Bulk Powder Enclosures

## Exclusive Features

- Retaining Ring provides a way to secure the source container liner bag.
- Retaining Band secures the source container liner bag to the retaining ring to seal the bag to the work surface.



# XPert® Bulk Powder Enclosures

## Exclusive Features



- Baffle pivots downward for easy cleaning.

# XPert® Bulk Powder Enclosures



- Digital airflow alarm for continual display of airflow velocity
- Hydraulic lift holds sash up while loading and unloading

# XPert® Bulk Powder Enclosures


## Access Opening



- Variable location
  - ✓ Right side
  - ✓ Center
  - ✓ Left side (not shown)

# XPert<sup>®</sup> Bulk Powder Enclosures

## Planning



- Understand procedures and equipment
- Arrange for minimal disruptions
- Consult your Safety Officer for protective equipment recommendations
- Transport containers to Bulk Powder Enclosure

# XPert® Bulk Powder Enclosures

## Start Up

- Turn on blower and purge system for several minutes before beginning procedure.
  - Check the baffle air slots for obstructions.



# XPert® Bulk Powder Enclosures

## Start Up



- Load with necessary equipment and supplies.
- Decontamination products may be useful.
- Ensure that sash is fully closed when working in enclosure.
- Position equipment so air flow is not obstructed.

# XPert® Bulk Powder Enclosures

## Loading



- Attach Safety Sleeve to source container and place on cart.
- Dock cart to Bulk Powder Transfer System.
- Push button to activate Lift System.
- Lift stops automatically when source container is at access height.



# XPert® Bulk Powder Enclosures

## Loading



- Pull Safety Sleeve through opening in work surface and secure on Retaining Ring with Retaining Band.
- Allow system to purge airborne contaminants from the work area.

# XPert® Bulk Powder Enclosures

## Loading



- Remove lid from source container.
  - ✓ Place inside enclosure, or
  - ✓ Decontaminate thoroughly and remove from enclosure
- Secure inner bags from source container to the Retaining Ring with Retaining Bands.

# XPert® Bulk Powder Enclosures

## Work Technique

- Keep all materials inside the lower air foil and perform all contaminated operations as far to the rear of the work area as possible.
- Segregate all clean and contaminated materials in the work area.
- Avoid using techniques or procedures that disrupt airflow patterns of the enclosure.



# XPert<sup>®</sup> Bulk Powder Enclosures

## Work Technique



- Keep all materials and hands inside enclosure until the end of the procedure and all materials have been wiped down.
- At end of wipe down, wipe down arm covers and gloves.
- Remove arm covers and gloves inside enclosure.
- Place waste materials in waste chute.

# XPert® Bulk Powder Enclosures

## Work Procedure



- Remove lid of source container and place to the side. Secure both outer and inner bags of source container to retaining ring with retaining band.
- Scoop powder into container.

# XPert® Bulk Powder Enclosures


## Work Procedure

When weighing procedure is completed, close container and place to the side. The user may either wipe container and lid after each weigh out or when all weigh outs are completed.



# XPert® Bulk Powder Enclosures

## Clean Up



- When procedure is completed, remove retaining band and place in disinfecting solution or discard in waste chute.
- Place inner bag in source container. Repeat for outer bag and place lid back on source container.
- Vacuum if necessary. With a decontamination solution such as IPE, wipe down source container lid, inside of safety sleeve, work surface and any other materials inside the enclosure.

# XPert® Bulk Powder Enclosures

## Clean Up



- Remove retaining band on safety sleeve.
- Pull safety sleeve up through opening in work surface and discard.
- With a decontamination solution, wipe down arm covers and gloves before removing and discarding in waste chute.



# XPert® Bulk Powder Enclosures

## Clean Up



- Lower source container until it comes to a stop.
- Unlatch Transfer Cart and return source container.
- Remove materials from inside enclosure and properly dispose of materials.
- Allow unit to purge for 2 minutes before turning off exhaust system.

# XPert® Balance Enclosures

## Accessories



- Two exhaust configurations
  - ✓ Top exhaust transition
  - ✓ Bottom exhaust transition
- Two ductwork options
  - ✓ 5" diameter flexible hose
  - ✓ 6" diameter PVC Duct

# XPert® Balance Enclosures

## Accessories



### Remote Blowers

- TEFC-type direct drive motor blower
- Corrosion-resistant phenolic-coated steel housing and impeller



- Belt drive motor blower
- Corrosion-resistant phenolic-coated steel housing and impeller
- Dry powder epoxy-coated base and weather cover

# XPert® Balance Enclosures

## Accessories

- Canopy Connection Kit
  - Remote Exhaust Fan
    - Ductwork
      - Dampers



# XPert® Balance Enclosures

## Accessories



- Kim wipe holder
- Spray bottle holder
- Ergonomic chair
- Ergonomic footrest



# XPert® Balance Enclosures


## Accessories



- 3-Shelf accessory kit
  - ✓ Reduces interior clutter
  - ✓ Shelves can rest on back wall
  - ✓ Does not affect airflow
- Printer shelf—inside or out
  - ✓ Utensil shelf
  - ✓ Utility shelf—disposables or calibration weight

# XPert® Balance Enclosures

## Accessories



- Trace Odor Carbon Filters for XPert System
- Carbon Filters
- Fluorescent Light
- HEPA Filters

# XPert® Balance Enclosures

## Accessories



- Solid Epoxy Work Surfaces
- Stainless Steel Work Surfaces





# XPert® Balance Enclosures

## Accessories



- Telescoping Base Stand
- Transfer Caddy
- Protector® Standard Storage Cabinets
- Telescoping Base Stands with Casters

# XPert® Balance Enclosures

## Conformity



- SEFA 1-2005
- ANSI Z9.5-1993
- Modified ASHRAE 110-1995
- UL 3101-1/61010-1 (115v models)
- CAN / CSA C22.2 No. 1010.1 (115v models)
- CE Conformity Marking (230v models)
- Safebridge Tested

# XPert® Balance Enclosures

## Additional Information

- Product literature
- Lunch-and-Learn sessions
- XPert video
- Trade shows
- 3-part specifications
- CAD drawings

