



NOTES:
 DUE TO DOOR DIMENSION WHERE THE WIDTH IS LARGER THAN THE HEIGHT,
 PASS THROUGH MUST BE MOUNTED SO THE DOORS HANG, LATCH & INTERLOCK
 PROPERLY. TERRAUNIVERSAL IS NOT RESPONSIBLE FOR
 IMPROPER FUNCTION FROM IMPROPER INSTALLATION

UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES TOLERANCES ON:		TITLE	
XX DECIMALS ± .03	X±	1"	PASS THROUGH, 36" W x 24" D x 24" H ACRYLIC
.XXX DECIMALS ± .010	X/X ±	1/64	
DRAWN: VEGA		DATE: 09/22/05	NEXT ASSY USED ON SCALE NONE PASS THROUGH SHEET 1 OF 1
SIGNATURES: TUCKER			
-1		127246	
CONFIG	DESCRIPTION	FINISH SPEC.	
MATERIAL SPEC.			
DATE: 04.26.05	APPRO: SV		
INITIAL RELEASE			
CHANGED RACKS			
MODIFIED HEIGHT			
DATE: 08.19.05	APPRO: SV		
DESCRIPTION			