



LTR	DESCRIPTION	DATE	APPD	COUNTS	DESCRIPTION
NC	INITIAL RELEASE	02.13.05	EDT	304SS 16 GA	FINISH SPEC. #4 BRUSH FINISH

UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES

TOLERANCES:	XX DECIMALS ±	X/2
	.XXX DECIMALS ±	X/X ±
DRAWN:	TUCKER	DATE
ENGR:	TUCKER	02.13.05

Terra Universal

PROFESSIONAL & COMMERCIAL
 A TERRA UNIVERSAL
 DISTRIBUTOR OF CORIAN PRODUCTS

TITLE: PASS THROUGH
 18"Wx18"Dx16"H
 SS, SDPVC WINDOW

SCALE: A

1992-11A

P/T SHEET 1 OF 1