



- NOTES:
1. SHOWN WITHOUT OPTIONAL MOUNTING BRACKETS. RECOMMENDED WALL CUTOUT DIMENSIONS ARE NOMINAL WIDTH AND HEIGHT DIMENSIONS PLUS 1/2"

LR#	DESCRIPTION	DATE	APPD	CONFIG	DESCRIPTION	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES TOLERANCES ARE:
NC	INITIAL RELEASE	09/13/04		-1		.XX DECIMALS ± .08 X±
					MATERIAL SPEC. 304 STAINLESS STEEL	.XXX DECIMALS ± X/X ±
					FINISH SPEC.	
DRAWN: GONZALEZ						DATE: 8/23/04
ENGR:						
						TITLE: PASS THROUGH 28.7" W x 34.5" D x 24" H Series 100 BioSafe EP SS
PROJECT: 2636-06A-2-58282 SHEET: 1 OF 2						