


NOTE-- THIS IS A CONTINUOUS-SEAM WELDED
304 STAINLESS STEEL DESSICCATOR CHAMBER

LIR	DESCRIPTION	DATE	APPD	CONTRIS	DESCRIPTION	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES
NG	INITIAL RELEASE	3/12/07	JY	304 SS	FINISH SPEC	TOLERANCES ON
				SPPVC		XX DECIMALS ± .03 X± 1"
						XXX DECIMALS ± .010 X/X ± 1/64
						DRAWN: YARKLE DATE: 3/12/07
						ENGINEER: _____
 <p>PROPRIETARY & CONFIDENTIAL TERRA UNIVERSAL DECLASSIFY OR CORRECT IF REQUIRED</p>						
<p>TYPE: PASS-THROUGH CHAMBER CONTINUOUS-SEAM; SS WALL MOUNT 24" X 12" X 36"</p>						
NEXT ASSY USED: 141220			SCALE: NONE		P/T: SHEET 1 OF 1	