


FRONT VIEW



RIGHT SIDE VIEW

NOTE - THIS IS A CONTINUOUS-SEAM WELDED
 316 STAINLESS STEEL PASS-THROUGH CHAMBER

L/R	DESCRIPTION	DATE	APPD	CONFIG	DESCRIPTION	FINISH SPEC
NG	INITIAL RELEASE			MATERIAL SPEC		
				316, SS		
UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES						
TOLERANCES ON						
XX	DECIMALS ±	.03	X±		1"	
XXX	DECIMALS ±	.010	X/X ±		1/64	
DRAWN:	SIGNATURES		DATE			
ENGR:	VEGA		12/07/05			
	TUCKER					
TITLE: CUSTOM PASS-THROUGH 316 SS, FLOOR MOUNT 48" X 48" X 60" H						
NEXT ASSY USED: 130823			SCALE: A		NONE	
C/R			SHEET 1		OF 1	