



ENGINEERING APPROVAL

ENGINEERING APPROVAL

CUSTOMER APPROVAL

LTR	DESCRIPTION	DATE	APPD	-1	CONFIG	DESCRIPTION
NC	INITIAL RELEASE	06.22.99	EDIT		MATERIAL SPEC.	FINISH SPEC.
					304SS BODY	
					304SS WINDOW	
UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES						
TOLERANCES ON:						
XX	DECIMALS ± .03	X ±	1"			
XXX	DECIMALS ± .010	X/X ±	1/64			
DRAWN:	TUCKER	DATE				
ENGR:	TUCKER	05.28.99				
		05.28.99				
PROPRIETARY & CONFIDENTIAL ALL RIGHTS RESERVED TERA UNIVERSAL DIVISION OF COMMERCE						
TITLE: PASS-THROUGH, SS FLOOR MOUNT, 36"Wx38"Dx60"H I.D.						
NEXT ASSY USED ON SCALE NONE PASSTHROUGH SHEET 1 OF 1						