



Invest Confidently with Metro 316 Stainless Steel Solutions.

Introducing Metro's 316 Stainless Steel Carts and Shelving

Engineered for applications demanding the highest standards, Metro's new Type 316 Stainless Steel products combine unparalleled durability, corrosion resistance, and cleanability. Whether you need mobile carts or stationary shelving, Metro provides the perfect solution.

Key Features of 316 Stainless Steel Carts and Shelving:

- **316 Stainless Steel Shelves & Corners:** Provide superior resistance to chlorides and harsh chemicals, ensuring durability in aggressive environments.
- **304 Stainless Steel Posts & Casters:** Durable construction with Metro's cart-washable polyurethane casters (2 swivel, 2 brake) for reliable mobility and ease of use.
- **Smart Design for Cleanability:** Flat, solid shelving with raised edges contains spills and minimizes contamination, while smooth surfaces allow for effortless cleaning.
- **Flexible and Versatile Configurations:** Adjustable shelving heights in 1" increments and modular design make it easy to adapt to changing needs, whether stationary or mobile.

Model No.	Size (in.)	(mm)	Shelf #	Approx. Pkd. Wt.
MCP316-1830FF-2	18x30	(457.2 x 762)	2	60.2 lbs. (27.3 kg)
MCP316-1830FF-3	18x30	(457.2 x 762)	3	66.7 lbs. (30.3 kg)
MCP316-2436FF-2	24x36	(609.6 x 914.4)	2	72.9 lbs. (33.1 kg)
MCP316-2436FF-3	24x36	(609.6 x 914.4)	3	79.1 lbs. (35.9 kg)

Load Capacity Tailored to Application:

- **Carts:** Load capacity may vary depending on casters and mobility factors but remains robust for high-demand applications.
- **Shelving Systems:** Shelves can support up to 800 lbs. (363 kg) for lengths 48" and shorter, with adjusted capacities for longer shelves (e.g., 600 lbs. for 60").



Applications for 316 Stainless Steel

Metro's 316 stainless steel products are ideal for industries requiring exceptional durability, cleanability, and resistance to harsh cleaning agents.



Laboratories & Cleanrooms

Handles aggressive cleaning agents and maintains sterility.



Healthcare

Perfect for operating rooms and sterile pharmacies.



Pharmaceuticals

Built for demanding, high-sanitization environments.

Model No.	Width/Length (in.)	(mm)	Approx. Pkd. Wt. (lbs.)	(kg)
1824NFS-T316	18x24	457x610	11	5.0
1830NFS-T316	18x30	457x760	14	6.4
1836NFS-T316	18x36	457x914	16	7.3
1842NFS-T316	18x42	457x1066	18	8.2
1848NFS-T316	18x48	457x1219	20	9.1
1860NFS-T316	18x60	457x1524	24	10.9
2424NFS-T316	24x24	610x610	15	6.8
2430NFS-T316	24x30	610x760	17	7.7
2436NFS-T316	24x36	610x914	19	8.6
2442NFS-T316	24x42	610x1066	21	9.5
2448NFS-T316	24x48	610x1219	24	10.9
2460NFS-T316	24x60	610x1524	29	14.0

Corrosive Resistance	Durability	Cleanability	Longevity
★★★★★	★★★★★	★★★★★	★★★★★

Discover the Metro Difference.

Uncompromising Quality

Type 316 stainless steel offers superior corrosion resistance, particularly against chlorides and aggressive chemicals. Its added molybdenum makes it ideal for demanding environments like cleanrooms, healthcare facilities, and laboratories that require frequent cleaning and sterilization. Products are designed for maximum durability and hygiene.

Long-Term Value

Reduced maintenance and extended product lifespan lower total cost of ownership.

Ready to elevate your standards?

Contact your Metro representative today to explore our Type 316 stainless steel carts and shelving. Together, we'll find the perfect solution for your laboratory or cleanroom needs.



We put space to work.

© 2025 InterMetro Industries Corporation, Wilkes-Barre, PA 18705

