

PRODUCT MODEL

Catalog # 5171810

8' Telescoping Base Stand with Leveling Feet and Shelf (non-welded), Optic White, with toe-locking casters

All Features & Specifications

Specifications

Applicable Products	Logic Biosafety Cabinets (september 2024 and Newer)
----------------------------	---

Works With These Products	CentriVap Pro Vacuum Concentrators, CentriVap Standard Vacuum Concentrators
----------------------------------	---

Built-in Options	Fixed Leveling Feet
-------------------------	---------------------