






ALL SHOPS



STATIC DISSIPATIVE
VIEWING WINDOW



MOUNTING BRACKET
(SHIP FIXED)

MOUNTING BRACKET
(SHIP LOOSE)

REVISIONS				MATERIAL SPEC.	FINISH SPEC.	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES		 Terra Universal Proprietary & Confidential All rights reserved by TERRA UNIVERSAL Disclosure or Copying Prohibited	TITLE	
RV#	DESCRIPTION	DATE	BY			TOLERANCES ON:			PASS-THROUGH: w/BRACKETS, SD PVC, WALL MOUNT 18"W x 18"D x 24"H	
NC	INITIAL RELEASE	1-19-17	JR			.XX DECIMALS ± .03 X" ± 1"				
○						.XXX DECIMALS ± .010 X/X ± 1/64				
						SIGNATURES			DATE	
						ISSUED BY: JUAN R			5/8/2019	
						REF #:			A	1992-40D
									SCALE	SHEET 1 OF 9