



# **CLEANROOM VINYL**

## **BENCH HEIGHT**

# 10 SERIES



AL10-VCR



PL10-VCR

Super Heavy Duty 450 lbs, Capacity, 3" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model    | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base        | Caster |
|----------|-----------------|-------------------|-----------------|-----------------|----------|-------------|--------|
| AL10-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 15" W x 12" H   | 20"      | Polished    | 1      |
| PL10-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 15" W x 12" H   | 20"      | Black Nylon | 1      |



AM10-VCR



PM10-VCR

Super Heavy Duty 450 lbs, Capacity, 3" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model    | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base              | Caster |
|----------|-----------------|-------------------|-----------------|-----------------|----------|-------------------|--------|
| AM10-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 16" W x 16" H   | 20"      | Polished Aluminum | 1      |
| PM10-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 16" W x 16" H   | 20"      | Black Nylon       | 1      |



AE10-VCR



PE10-VCR

Super Heavy Duty 450 lbs, Capacity, 3" High Resilient Density Foam for Better Comfort, Rubber Edge  
Banding Protection on Seat & Back

| Model    | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base              | Caster |
|----------|-----------------|-------------------|-----------------|-----------------|----------|-------------------|--------|
| AE10-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 19" W x 16.5" H | 20"      | Polished Aluminum | 1      |
| PE10-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 19" W x 16.5" H | 20"      | Black Nylon       | 1      |

# 20S SERIES



AM20S-VCR



PM20S-VCR

Super Heavy Duty 450 lbs, Capacity, 2.5" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model     | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base              | Caster |
|-----------|-----------------|-------------------|-----------------|-----------------|----------|-------------------|--------|
| AM20S-VCR | Cleanroom       | 21"-31"           | 18" W x 18" D   | 16" W x 16" H   | 18"      | Polished Aluminum | 1      |
| PM20S-VCR | Cleanroom       | 21"-31"           | 18" W x 18" D   | 16" W x 16" H   | 18"      | Black Nylon       | 1      |



AL20S-VCR



PL20S-VCR

Super Heavy Duty 450 lbs, Capacity, 2.5" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model     | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base              | Caster |
|-----------|-----------------|-------------------|-----------------|-----------------|----------|-------------------|--------|
| AL20S-VCR | Cleanroom       | 21"-31"           | 18" W x 18" D   | 15" W x 12" H   | 18"      | Polished Aluminum | 1      |
| PL20S-VCR | Cleanroom       | 21"-31"           | 18" W x 18" D   | 15" W x 12" H   | 18"      | Black Nylon       | 1      |

# 20M SERIES



AM20M-VCR



PM20M-VCR

Super Heavy Duty 450 lbs, Capacity, 3" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model     | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base        | Caster |
|-----------|-----------------|-------------------|-----------------|-----------------|----------|-------------|--------|
| AM20M-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 16" W x 16" H   | 20"      | Polished    | 1      |
| PM20M-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 16" W x 16" H   | 20"      | Black Nylon | 1      |



AL20M-VCR



PL20M-VCR

Super Heavy Duty 450 lbs, Capacity, 3" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model     | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base        | Caster |
|-----------|-----------------|-------------------|-----------------|-----------------|----------|-------------|--------|
| AL20M-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 15" W x 12" H   | 20"      | Polished    | 1      |
| PL20M-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 15" W x 12" H   | 20"      | Black Nylon | 1      |



AE20M-VCR



PE20M-VCR

Super Heavy Duty 450 lbs, Capacity, 3" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model     | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base              | Caster |
|-----------|-----------------|-------------------|-----------------|-----------------|----------|-------------------|--------|
| AE20M-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 19" W x 16.5" H | 20"      | Polished Aluminum | 1      |
| PE20M-VCR | Cleanroom Vinyl | 21"-31"           | 20" W x 18" D   | 19" W x 16.5" H | 20"      | Black Nylon       | 1      |

# 20W SERIES



AE20W-VCR



PE20W-VCR

Super Heavy Duty 450 lbs, Capacity, 3" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model     | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base              | Caster |
|-----------|-----------------|-------------------|-----------------|-----------------|----------|-------------------|--------|
| AE20W-VCR | Cleanroom Vinyl | 21"-31"           | 21.5" W x 18" D | 19" W x 16.5" H | 20"      | Polished Aluminum | 1      |
| PE20W-VCR | Cleanroom Vinyl | 21"-31"           | 21.5" W x 18" D | 16" W x 16.5" H | 20"      | Black Nylon       | 1      |



AL20W-VCR



PL20W-VCR

Super Heavy Duty 450 lbs, Capacity, 3" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model     | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base              | Caster |
|-----------|-----------------|-------------------|-----------------|-----------------|----------|-------------------|--------|
| AL20W-VCR | Cleanroom Vinyl | 21"-31"           | 21.5" W x 18" D | 16" W x 16" H   | 20"      | Polished Aluminum | 1      |
| PL20W-VCR | Cleanroom Vinyl | 21"-31"           | 21.5" W x 18" D | 16" W x 16" H   | 20"      | Black Nylon       | 1      |



AM20W-VCR



PM20W-VCR

Super Heavy Duty 450 lbs, Capacity, 3" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model     | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base              | Caster |
|-----------|-----------------|-------------------|-----------------|-----------------|----------|-------------------|--------|
| AM20W-VCR | Cleanroom Vinyl | 21"-31"           | 21.5" W x 18" D | 16" W x 16" H   | 20"      | Polished Aluminum | 1      |
| PM20W-VCR | Cleanroom Vinyl | 21"-31"           | 21.5" W x 18" D | 16" W x 16" H   | 20"      | Black Nylon       | 1      |

## 45 SERIES



A45-VCR



P45-VCR

Super Heavy Duty 300 lbs, Capacity, 2.5" High Resilient Density Foam for Better Comfort, Rubber Edge Banding on Seat and Back for a good Protection of Fabric and Vinyl.

| Model   | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base              | Caster |
|---------|-----------------|-------------------|-----------------|-----------------|----------|-------------------|--------|
| A45-VCR | Cleanroom Vinyl | 21"-31"           | 18"W x 17"D     | 15.5"W x 13"H   | 18"      | Polished Aluminum | 1      |
| P45-VCR | Cleanroom Vinyl | 21"-31"           | 18"W x 17"D     | 15.5"W x 13"H   | 18"      | Black Nylon       | 1      |

## 50 SERIES



50-VCR

Super Heavy Duty 300 lbs, Capacity, 2.5" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & Back

| Model  | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base        | Caster |
|--------|-----------------|-------------------|-----------------|-----------------|----------|-------------|--------|
| 50-VCR | Cleanroom Vinyl | 19"-27"           | 18" W x 17" D   | 15" W x 12" H   | 18"      | Black Nylon | 1      |

# 60 SERIES



A60-VCR



P60-VCR

Super Heavy Duty 300 lbs, Capacity, 2.5" High Resilient Density Foam for Better Comfort, Rubber Edge Banding Protection on Seat & back

| Model   | Upholstery Type | Height Adjustment | Seat Dimensions | Back Dimensions | Footring | Base        | Caster |
|---------|-----------------|-------------------|-----------------|-----------------|----------|-------------|--------|
| A60-VCR | Cleanroom Vinyl | 21"-31"           | 18"W x 17"D     | 15.5"W x 13"H   | 18"      | Polished    | 1      |
| P60-VCR | Cleanroom Vinyl | 21"-31"           | 18"W x 17"D     | 15.5"W x 13"H   | 18"      | Black Nylon | 1      |

Industrial Seating™ offers clean room seating designed to maintain reduced particulate levels as specified in federal standards 209E, Class 10/.2 micron and or class 100/.5 micron **Sealed System:** On Industrial Seating™ cleanroom seating, the air passed through three 99.99% effective .2 micron filter systems. These filters allow air within the sealed cushions to be expelled under minimal pressure. This filtering system first "strains out" larger particles, then traps all particles down to below .02 micron size. Permanent adhesives are used to doubled seal edges of all cover materials used on Industrial Seating™ class 10 and class 100 clean room seating.

## PreFixx Group II

- 100 Vinyl Face, Poly/cotton Backing, • Protective Coating for Increased Stain Resistance, • Abrasion 400,000 Double Rubs
- Meets or Exceeds California Tech Bulletin 117, Boston Fire Dept. IX-1, • Federal Standard A-A-2950, • UFAC Fabric Class 1, • FMVSS-302

|          |          |              |           |          |          |          |
|----------|----------|--------------|-----------|----------|----------|----------|
|          |          |              |           |          |          |          |
| Blue-211 | Navy-221 | Burgundy-261 | Black-251 | Grey-231 | Grey-232 | Grey-233 |

# OPTIONS



DLC



ST



FLC

| Controls | Description  |
|----------|--|
| DLC      | Height and Back Tilt (This control has two levers for seat height and back angle adjustment and back height adjustment with knob)                                  |
| ST       | Height and Seat and Back Tilt ( This control has two levers, one for height adjustment and the other for Seat and Back tilt adjustment and back height with knob)  |
| FLC      | Three Lever for Individual Lift ( Height lever adjustment, seat tilt adjustment back tilt adjustment and height back adjustment with knob and a seat tension knob) |



FR04



FR05



FR06

| Footring | Description   |
|----------|---|
| FR04     | Knobless Nylon Footring 18"<br>The Footring is Loosened and Tightened by Turning the Outer Ring. (option for models 20S, 45, 50 and 60) |
| FR05     | Knobless Aluminum with Steel Ring 20"<br>The Footring is Loosened and Tightened by Turning the Outer Footring. (option                  |
| FR06     | Nylon Center with Steel Ring, Footring 18" (option for models 20S, 45, 50 and 60)   |



TF



TA



TAA

| Arms | Description                              |
|------|--|
| TF   | Fixed Height and Width, Black            |
| TA   | Fixed Height and Adjustable Width, Black |
| TAA  | Adjustable Height and Adjustable Width   |



1B



1U



1UB



2

| Caster | Description  |
|--------|--|
| 1B     | Braking Caster , Recommended where Chairs may roll to easily when is Not Occupied  |
| 1U     | Soft Urethane Caster, Recommended on Tile, Cement and other Hard Surfaces          |
| 1UB    | Soft Urethane Braking Caster, Recommended for reduce Mobility when is Not Occupied |
| 2      | Nylon Glide, Recommended for use in areas where requiring limited mobility         |



TF-BX-C



TA-BX-C



TAA-BX-C

| Arms     | Description   |
|----------|---|
| TF-BX-C  | Fixed Height and Width, with Chrome Arm Frames & ESD Arm Pads                 |
| TA-BX-C  | Fixed Height and Adjustable Width, with Chrome Arm Frames & ESD Arm Pads      |
| TAA-BX-C | Adjustable Height and Adjustable Width, with Chrome Arm Frames & ESD Arm Pads |



9B



9U



9UB



11

| Caster | Description  |
|--------|--|
| 9B     | ESD Braking Caster , Recommended where Chairs may roll to easily when is       |
| 9U     | ESD Soft Urethane Caster, Recommended on Tile, Cement and other Hard           |
| 9UB    | ESD Soft Urethane Braking Caster, Recommended for reduce Mobility when is      |
| 11     | ESD Nylon Glide, Recommended for use in areas where requiring limited mobility |