







SIGN OFF



Fire-rated door for 90-minute fire exposure



REVISIONS			
RV#	DESCRIPTION	DATE	BY
NC	INITIAL RELEASE	8/4/22	MM

MATERIAL SPEC.
316L SS

FINISH SPEC.

UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES	
TOLERANCES ON:	
.XX DECIMALS ± .03	X' ± 1"
.XXX DECIMALS ± .010	X/X ± 1/64
SIGNATURES	DATE
ISSUED BY: MARCO MELENDEZ	8/5/2022
REF #:	-



PROPRIETARY & CONFIDENTIAL
ALL RIGHTS RESERVED BY
TERRA UNIVERSAL
DISCLOSURE OR COPYING PROHIBITED

TITLE
PASSTHROUGH; SMART, FIRE RATED DOOR CLEANMOUNT CLEANSEAM 24"W X 24"D X 24"H ID 316L SS WALL MOUNT
A 2635-10B-2-316-FR
SCALE
SHEET 1 OF 1