



23.500
OD CHAMBER



HINGE


SS LIFT LATCH

SS CATCH

MOUNT BRACKET

COVER MOUNT
BRACKET


ACCESSE PANEL



17.589

18.000

28.560

WALL CUT OUT
W 27.308 X H 18.385ID 18.000
CHAMBER

20.653

TEMPERED GLASS
PA07456
(X2)

FLUSH MOUNT

NUT

HANDLE
PA06513-16
(X2)

GUIDE

SCREW

16.000

12.200

17.560

REVISIONS				MATERIAL SPEC.	FINISH SPEC.	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES				Terra Universal PROPRIETARY & CONFIDENTIAL ALL RIGHTS RESERVED BY TERRA UNIVERSAL DISCLOSURE OR COPYING PROHIBITED	TITLE	
RV#	DESCRIPTION	DATE	BY			TOLERANCES ON:					A SCALE	2636-76C SHEET 1 OF 20
NC	INITIAL RELEASE	3-1-17	JR			.XX DECIMALS ±.03 X" ± 1"						
						.XXX DECIMALS ±.010 X/X ±1/64						
						SIGNATURES DATE						
						ISSUED BY: JUAN R 3-1-17						
						APPROVED BY:						
				REF #:								

SCALE

SHEET 1 OF 20