


Specifications	Model Number				Representative Model (tested)
	TSX505GA				
	Recommended cabinet temperature setpoint		Capacity m ³		
	5°C		0.16 (5.5 cubic feet)		
	Application				
Vaccine storage					
Program	Temperature requirements				
	Steady state	Short door openings	Long (3 min.) door opening recovery		
NSF 456	2°C < T < 8°C	2°C < T < 8°C	in < 15min, < 8°C		TSX505GA

Empty, TSX505GA, Under-counter Refrigerator @ 5°C, NSF test sequence @ 22°C



Loaded, TSX505GA, Under-counter Refrigerator @ 5°C, NSF test sequence @ 22°C

