



REVISIONS			
RV#	DESCRIPTION	DATE	BY
NC	INITIAL RELEASE	01/24/18	

MATERIAL SPEC.	304 SS
FINISH SPEC.	1/2" TEMPERED GLASS

UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES			
TOLERANCES ON:			
.XX DECIMALS ±.03	X' ± 1"		
.XXX DECIMALS ±.010	X/X' ± 1/64		
SIGNATURES		DATE	
ISSUED BY:		01/24/18	
APPROVED BY:			
REF #:			



Terra Universal

PROPRIETARY & CONFIDENTIAL
ALL RIGHTS RESERVED BY
TERRA UNIVERSAL
DISCLOSURE OR COPYING PROHIBITED

TITLE	Window, Framed; Single-Sided, Tempered Glass, BioSafe Overall Dimensions 36"W x 48"H Wall Thickness 5"
SCALE	6611-00-1204-116181
SHEET	1 OF 1