

Biohazard Bin XL

Lab Waste Container

- 23 gallon capacity
- Foot pedal operated lid
- Printed biohazard symbol
- Bag holder support notches



This full-size 23 gallon (87L) biohazard bin is perfect for all of the infectious or contaminated waste produced in medical facilities and laboratories on a daily basis. Permanently marked with the international biohazard symbol, there is no need for additional stickers or labels. The foot-pedal operated lid allows for hand free operation and reduces contamination risks. Notch hooks on the outside of the basket rim make installing bags simple and easy.

Specifications

Color:	Red - Biohazard
Dimensions (in.):	19 x 19 x 33in.*
Dimensions (cm):	48 x 48 x 84cm*
Capacity:	23 gallons / 87L

* including lid

For Lab & Medical Use

Cat. No	Description	Qty	List Price
A8002B	Biohazard Bin XL, large, 23 gallon, with foot pedal operated lid, 19 x 19 x 33in. (LxWxH)	1	\$189.00

Other available sizes...

Cat. No	Description	Qty	List Price
A8000B	Biohazard Bin, medium, 13 gallon, with foot pedal operated lid, 10 x 14 x 19in. (LxWxH)	1	\$74.79
A8401B	BenchBin™ benchtop biohazard bin, 7 x 5 x 7in. (LxWxH)	1	\$39.80