



Website Prices Honored Worldwide

Critical Environment Solutions

Holder; Safety Glasses, 304 SS, 10.1"W x 6.25"D x 26"H, 24 Compartments, Wall Mount



Price: **\$1,140.00 EA**

Part #: 4949-04

Usually ships in 1-3 days

Manufacturer: Terra Universal, Inc.

U.S. orders include packaging.

Product URL: <https://www.terrauniversal.com/product/4949-04>

24-chamber stainless steel safety glasses holder with swing-out hinges.

Save space and organize safety glasses

- Includes keyhole cut-out for wall mounting; due to variations in wall construction and mounting codes, fasteners are not included.
- Hinged front compartment can be opened with draw latch for rear access to load glasses
- Holder can be lifted off the wall for thorough cleaning; drop hinge pins back into place to reattach
- Stainless steel is resistant to most alcohol-based cleaners and common lab chemicals
- BioSafe design includes smooth, easily cleanable edges

Package Type: Box

Packaged Weight: 39 lbs

Product Dimensions, in. 10.1 W, 6.25 D, 26 H

Package Dimensions, in. 13 W, 9 D, 33 H

The price and ship time shown may have changed since this document was printed or saved.
Current product information can be found on the [TerraUniversal.com](https://www.terrauniversal.com) website or by calling Terra at 714-578-6000.