

TERRA UNIVERSAL Critical Environment Solutions

Website Prices Honored Worldwide

Stainless Steel Flat Top Table; 60" W x 30" D x 34.5" H, Sliding Door



Price: \$2,011.00 EA

Part #: 5656-11A

Usually ships in 7-13 days

Manufacturer:

U.S. orders include packaging.

Product URL: <https://www.terrauniversal.com/product/5656-11A>

All doors and panels have safety edges and are free from bacteria-gathering crevices

- Material: 14 gauge type 304 stainless steel
- Contaminant-free and easy to maintain
- Flat top configuration
- All-welded design
- Stainless steel adjustable bullet feet
- Rear edges are turned down on 90° angle for alignment with walls or other equipment

Electrical: 12 V V, A, W, P, Hz

Package Type: Crate

Packaged Weight: 287 lbs

Product Dimensions, in. 60 W, 30 D, 34.5 H

Package Dimensions, in. 65 W, 35 D, 39 H

The price and ship time shown may have changed since this document was printed or saved.
Current product information can be found on the [TerraUniversal.com](https://www.terrauniversal.com) website or by calling Terra at 714-578-6000.